

UDI Laser Marking Services

Laser marking services are now available at our
National Repair Centre.

Quantity discounts may apply.

On-Site (mobile lab) or **In-Site** (at your facility)
laser marking services may be available for
larger projects. **Contact us for more details!**

575 Wentworth Street East, Unit 29
Oshawa, Ontario L1H 3V8
Tel: 905.721.0740 / 1.800.387.2281
Email: info@mobileinstrument.ca
www.mobileinstrument.ca

**One Year
Readability
Guarantee!**

The use of **Unique Device Identifiers** (UDIs) in conjunction with instrument tracking systems is quickly becoming an industry standard for medical devices.

MISR has a solution to mark a durable DataMatrix on a wide variety of materials including stainless steel, titanium, anodized titanium, aluminum, plastics, and rubber using the innovative **Nux-20 Laser Marker** from NuTrace.

**Call us today at 1.800.387.2281
for a free quotation and a
complete list of our capabilities!**

MISR laser marks are guaranteed to read for one year or they will be re-applied at no charge. MISR will remove at no charge a defective laser mark from any item where removal is possible (e.g. stainless steel); however, MISR assumes no responsibility for further repair or replacement of any device where a laser mark is defective.

575 Wentworth Street East, Unit 29
Oshawa, Ontario L1H 3V8
Tel: 905.721.0740 / 1.800.387.2281
Email: info@mobileinstrument.ca
www.mobileinstrument.ca