


Comprehensive care for your surgical instruments

INSTRUMENT COLOUR CODING



Colour code to personalize your instruments for easy identification, grouping and locating of specific instruments or special instrument sets. ID by coating one ring, two rings, or the proximal end of thumb forceps.

The advantages of MISR Colour Coding:

- Does not interfere with the operation of the instrument
- Long lasting - will withstand steam sterilization
- Superior to standard taping - tape will start to crack and chip off over time
- Duralon "J" Type 11 Nylon Powder - Medical Grade
- Colours available: black, blue, yellow, green and red

For more information, please contact your On-Site MISR Repair Technician

