


## Flexible Endoscope Repairs

Mobile Instrument Service & Repair  
now offers repair services for ALL  
Flexible Endoscopes



For more information, contact:  
**Taylor Bailie**, Business Development  
905.244.1042 | [taylorbailie@mobileinstrument.ca](mailto:taylorbailie@mobileinstrument.ca)

29-575 Wentworth Street East  
Oshawa, Ontario L1H 3V8  
Tel: 905.721.0740 / 1.800.387.2281  
Email: [info@mobileinstrument.ca](mailto:info@mobileinstrument.ca)  
[www.mobileinstrument.ca](http://www.mobileinstrument.ca)



## **Major & Minor Repairs Available:**

- Air / water clog repairs and nozzle replacement
- Bending sheath replacement
- Angulation adjustment
- Biopsy channel
- Fluid invasion
- Suction cylinder & channel replacement
- Light bundle replacement
- Light guide umbilical replacement
- Insertion tube replacement
- Image bundle replacement

## **All Repairs Include:**

- Thorough diagnostic evaluation of damage
- Leak & shock testing
- Disassembling & properly replacing all components
- The best replacement components available
- Stringent inspection
- Fast turnaround time
- 6 month warranty
- Flat rate & tiered pricing