

INSTRUMENT COLOUR CODING

Mobile Instrument Service & Repair
offers **Instrument Colour Coding**

ADVANTAGES OF MISR COLOUR CODING:

- Does not interfere with the operation of the instrument
- Long lasting - will withstand steam sterilization
- Superior to standard taping - tape will start to crack and chip off over time
- Duralon "J" Type 11 Nylon Powder - Medical Grade
- Colours available: black, blue, yellow, green and red



For more information, contact us at:
905.721.0740 | 800.387.2281
info@mobileinstrument.ca



575 Wentworth Street East, Unit 29, Oshawa, ON L1H 3V8
Tel: 905.721.0740 | 800.387.2281
info@mobileinstrument.ca | www.mobileinstrument.ca