








Comprehensive care for your surgical instruments









## FIBER OPTIC CABLES

Fiber optic cables are available in standard and stainless steel reinforced Monocoil. Insert the appropriate code under each category to create the product code for your order. **For a Monocoil/reinforced cable add "XX"** to the end of your product code below.

	A Light Source Fitting	B Bundle Size	C Instrument Fitting	D Cable Length	Optional Monocoil Cable	Monocoil
FOC-					XX	 <p style="font-size: 8px;">Add monocoil underneath the sheath to enhance durability and protect fiberoptic stands.</p>

Code	A Light Source Fitting	
1	ACMI 	
2	WOLF Dyonics 	
3	Wolf Extended 	
4	Karl Storz 	
5	Olympus 	
6	Pilling 	

Code	B Bundle Size
1	2mm
4	3.5mm
7	5mm

Code	C Instrument Fitting Styles	
1	ACMI Female 	
2	ACMI Female Clip On 	
3	ACMI Male 	
4	Pilling 	
5	Karl Storz 	
6	Wolf Dyonics 	
7	Wolf Clip On 	
8	Olympus 	

Code	D Cable Length
4	7.5 ft.
7	10 ft.

Standard Colour: Gray

Custom length, size and colour available upon request.

For more information, please contact your On-Site MISR Repair Technician  
[info@mobileinstrument.ca](mailto:info@mobileinstrument.ca)    1.800.387.2281

